



Alcohol Burners



Model 3040

- Glass Lamp w/ Nickel-Plated Brass Cap
- Large 150ml Capacity
- Ideal For Students

VEE GEE Model 3040 Alcohol Burners feature a sturdy, heavy-walled glass design with a large capacity of 150ml. Fitted with a nickel-plated brass cap used to extinguish the flame and to prevent the evaporation of alcohol. Either methyl or ethyl alcohol can be used with the lamp.

Specifications:

Capacity: 150ml
Diameter: 80mm
Height: 100mm
Weight: 225g

Wick. Long-burning, easily replaceable.

Wick Holder. Durable, nickel-plated brass.

Burner. Large 150ml capacity; heavy-walled glass.

Cap. Durable, nickel-plated brass.



Ordering Information

Cat. No.	Description	Capacity	Diameter	Height	Weight
3040	Alcohol Burner, Glass	150ml	80mm	100mm	240g
3040-W	Replacement Wicks (pk/12)				



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veege.com • e-mail: sales@veege.net