



S10 D - 7G - KS - 110 Dispersing element

Technical Data

Volume range min. (H2O) [l]	0.001
Volume range max. (H2O) [l]	0.040
Stator diameter [mm]	7
Rotor diameter [mm]	4.8
Circumferential speed max. [m/s]	7.5
Immersion depth min. [mm]	20
Immersion depth max. [mm]	90
Shaft length [mm]	110
Material in contact with media	Polycarbonat (PC), Polysulfon (PSU)
Working temperature max. [°C]	100
Sterilization methods	autoclavable
Ident. No.	3433312