



S 25 D - 10 G - KS Dispersing element (10 pcs)

incl. 1 Disposable tube
 Plastic materials used approved by FDA.

Technical Data	
Volume range min. (H2O) [l]	0.01
Volume range max. (H2O) [l]	0.1
Stator diameter [mm]	10
Rotor diameter [mm]	6.75
Circumferential speed max. [m/s]	8.5
Immersion depth min. [mm]	15
Immersion depth max. [mm]	85
Shaft length [mm]	150
Material in contact with media	Polycarbonat (PC), Polyetheretherketon (PEEK)
Working temperature max. [°C]	100
Sterilization methods	autoclavable
Ident. No.	3452200