



### S 25 D - 14 G - KS Dispersing element (10 pcs)

incl. 1 Disposable tube  
 Plastic materials used approved by FDA.

Technical Data	
Volume range min. (H2O) [l]	0.01
Volume range max. (H2O) [l]	0.5
Stator diameter [mm]	14
Rotor diameter [mm]	9.5
Circumferential speed max. [m/s]	12
Immersion depth min. [mm]	15
Immersion depth max. [mm]	85
Shaft length [mm]	150
Material in contact with media	Polycarbonat (PC), Polyetheretherketon (PEEK)
Working temperature max. [°C]	100
Sterilization methods	autoclavable
<b>Ident. No.</b>	<b>3452100</b>